

Analysis of X-Ray Nike Spectra from Highly-Charged High-Z Ions

Y. Aglitskiy¹, J.L. Weaver², M. Karasik², V. Serlin², S.P. Obenshain², and Yu. Ralchenko³

¹Leidos, Inc., 11951 Freedom Drive, Reston, VA 20190, USA

²Naval Research Laboratory, 4555 Overlook Ave. SW, Washington, DC 20375, USA

³National Institute of Standards and Technology, Gaithersburg, MD 20899, USA

Joint ICTP-IAEA Advanced School and Workshop on Modern Methods in Plasma Spectroscopy, Trieste, Italy, 16-27 March 2015

Supported by DOE NNSA



Nike Laser Target Facility




Current Nike operating parameters:
 248 nm central wavelength, KrF laser
 Echelon-free ISI beam smoothing, 1-2 THz bandwidth
 Total laser energy: 2-3 kJ
 44 main beams; 12 backlighter beams
 Main beam focal profile: 0.2-1.2 mm FWHM
 Main beam pulse duration: 0.4-12 ns
 Main beam overlapped intensity:




typical $\sim 1 \times 10^{14}$ W/cm² (4 ns pulse, 0.75 mm FWHM spot)
 up to $1-2 \times 10^{15}$ W/cm² (0.4 ns pulse, 0.2 mm FWHM spot)

Sample Target Focal Profile



Sample Spike Pulse



NRL Nike *spectroscopy* program

- Successful development and application of state-of-the-art monochromatic x-ray imaging and spectroscopic techniques for study of IC plasma instabilities
 - Analysis of Richtmeyer-Meshkov and Raleigh-Taylor instabilities
 - X-ray imaging with spherically bent crystals
 - Measurement of relativistic electron transport
 - Use of bent Bragg-Fresnel lenses for imaging




Goals

- Test spectroscopic equipment for NIF
 - Convex crystal spectrometer
 - Spherical imaging spectrometer
- Identify **new spectral lines** for plasma diagnostics
- Provide **accurate x-ray spectra** for validation of NLTE plasma kinetics codes (cf. NLTE Code Comparison Workshops)




Generic outline of convex crystal spectrometer



Crystal : mica
 $2d = 20 \text{ \AA}$
 Working orders of reflection:
 $n = 1, 2, 3, 5, 7, \dots$
 $\lambda/\Delta\lambda > 1000$




Nike half-moon convex crystal spectrometer




Crystal : mica
 $2d = 20 \text{ \AA}$
 Working orders of reflection:
 $n = 1, 2, 3, 5, 7, \dots$
 $\lambda/\Delta\lambda > 1000$




Nike half-moon convex "full range" crystal spectrometer with the reference spectra of H- and He-like Al and Si and Ne-like Y, Zr, Nb, Mo and Ag




Spherical crystal spectrometer



W spectra: **irradiated foil (few tens of μm)** : 7-8 shots



Measured elements (as of Sep 2014)




Primary observed ions

- **H-like (1s) and He-like (1s²) ions**
 - Mg, Al, Si
- **Ne-like (2s²2p⁶) and adjacent ions**
 - Y, Mo, Ag,...
- **Ni-like (3d¹⁰) and adjacent ions**
 - W, Pt, Au,...

Element	Seq.	IP (eV)
Si	H	2673
Mo	Ne	4259
Ag	Ne	5558
Hf	Ni	3741
W	Ni	4057
Pt	Ni	4715
Au	Ni	4888



Cu and Si spectra




Collisional-radiative modeling

- **Flexible Atomic Code** (Gu, 2003)
 - Relativistic model potential w/ some QED
 - Level energies
 - Radiative and autoionization probabilities
 - Collisional cross sections
 - Cross sections fit and databased
- **NOMAD** (Ralchenko & Maron, 2001)
 - Time-dependent collisional-radiative code for non-Maxwellian plasmas
 - CX, laser pumping, parabolic states,.....
 - Steady state
 - Opacity (escape factor for populations, radiative transfer for spectra)
 - IP lowering (Stewart & Pyatt)
 - Typically [Ga]-[Fe] ions, ~35,000 levels


Energy level structure for Ni-like ions




Isoelectronic sequence of Ni-like ions



Au shot, Nov 2013




Au shot, Nov 2013




NIST
National Institute of
Standards and Technology
U.S. Department of Commerce

Au shot, Nov 2013




NIST
National Institute of
Standards and Technology
U.S. Department of Commerce

Au shot, Nov 2013




NIST
National Institute of
Standards and Technology
U.S. Department of Commerce

Au: convex Nov 2013 vs focusing May 2014




NIST
National Institute of
Standards and Technology
U.S. Department of Commerce

Mo reference spectrum




NIST
National Institute of
Standards and Technology
U.S. Department of Commerce

Mo reference spectrum




NIST
National Institute of
Standards and Technology
U.S. Department of Commerce

Au: focusing spectrometer : one shot




NIST
National Institute of
Standards and Technology
U.S. Department of Commerce

Au: focusing spectrometer, May 2014




NIST
National Institute of
Standards and Technology
U.S. Department of Commerce

Au: focusing spectrometer, May 2014




$T_e = 1900 \text{ eV}$
 $n_e = 10^{21} \text{ cm}^{-3}$
 $R = 0.2 \text{ mm}$

Au: focusing spectrometer, May 2014




$T_e = 1900 \text{ eV}$
 $n_e = 10^{21} \text{ cm}^{-3}$
 $R = 0.2 \text{ mm}$

Au: focusing spectrometer, May 2014




$T_e = 1900 \text{ eV}$
 $n_e = 10^{21} \text{ cm}^{-3}$
 $R = 0.2 \text{ mm}$

Au: focusing spectrometer, May 2014




$T_e = 1900 \text{ eV}$
 $n_e = 10^{21} \text{ cm}^{-3}$
 $R = 0.2 \text{ mm}$


Au: focusing spectrometer, May 2014



$T_e = 1900 \text{ eV}$
 $n_e = 10^{21} \text{ cm}^{-3}$
 $R = 0.2 \text{ mm}$




Pt : focusing, May 2014



E2 lines:
 $3p^6 3d^{10} - 3p^5 3d^{10} 4f$
 $3p^6 3d^{10} 4s - 3p^5 3d^{10} 4s 4f$

Ionization distribution: NOMAD vs FLYCHK



<http://nlte.nist.gov/FLY>

Conclusions

- New spectroscopic system on Nike KrF laser allows high-accuracy measurements of x-ray spectra for low-, mid- and high-Z elements (Mg to Bi)
- Typical electron temperatures and densities for irradiated foil experiments are about 2000 eV and 10^{21} cm^{-3}
- The measured spectra make up an accurate testbed for NLTE kinetic codes (NLTE-9?)
- Collaborations are **very** welcome

